


ELISA Plates

- C bottom & F bottom are both available.
- SBS footprint could work with many brands of automatic working station.
- Detachable plate is convenient for researchers to transfer sample for different treatment.





ORDER INFORMATION

Size	Туре	Material	Packaging	Cat.#	Quantity per Case	Size/Case(CM)			Weight/
						L	W	н	Case(Kg)
$0.2mL \times 96$, C bottom, High binding	Detachable	ABS+PS	10PCS/Bag	EP00001H	10 Bags/Case	50	31	16	6.00
0.2mL x 96, C bottom, Middle binding	Detachable	ABS+PS	10PCS/Bag	EP00001M	10 Bags/Case	50	31	16	6.00
0.2mL x 96, F bottom, High binding	Detachable	ABS+PS	10PCS/Bag	EP00006H	10 Bags/Case	50	31	16	6.00